

VORTEX MIXER (NB-105V)



Optional Accessory

NB-105V01	Mixing Cup
NB-105V02	Platform Head
NB-105V03	Micro plate Kit
NB-105V04	Multi sample plate Kit



Features

- Automatic start when vortex is touched on the mixing cup.
- Constant on after starting
- Adjustable RPM up to 3000rpm
- Wide option for accessories
- High quality design for safety
- CE marked

Specification

Items	Unit	NB-105V
Speed range	rpm	0 to 3000rpm
Operation Mode		Touch On & Constant On
Material		Aluminum Die-casting & ABS
Head	mm	Rubber Mixing Cup
Dimensions	kg	130(W)x155(D)x155(H)mm 3Kg
Power	V/Hz	AC220V, 0.14A, 60Hz